


ho CryoPen[®] o

Beauty Renewed



CryoPen[®] o

OPTIMAL - All-purpose, powerful cryosurgery device for regular use.



- 1-3mm applicator
- 2-6mm applicator
- 4-8mm applicator

Delivered with: 1x CryoPen O pen • 1x 16g N₂O cartridge • 3 applicators. Order additional 8g or 16g cartridges separately.

BETTER TREATMENTS

Effective, no-anesthetic removal of small or large lesions, in just a few seconds.

FASTER RECOVERY

Perfect results for complete removal of all benign skin lesions.

EXTREME PRECISION

Microfine jet for precise control.
Safe for use on delicate skin or near the eyes.

EXTREME POWER

N₂O at high pressure (55 bars) and ultralow temperature (-89°C) for effective cell destruction.

BETTER BUSINESS

Efficient, practical and extremely cost effective. For a quick return on investment. Ideal for medical professionals and estheticians.