

"Dual Pack" Unshielded Trocar & Cannula (Single Patient Use)



The Femcare-Nikomed "Dual Pack" Unshielded 5mm Trocar and Cannula is the perfect secondary port and includes a leak free seal which maintains pneumoperitoneum throughout surgery.

- Bladeless tip for easy insertion
- Silicon seal for maintaining pneumoperitoneum
- Cost effective

Product Description	Product Code	Pack Size
Unshielded Trocar & Cannula Dual Pack: 5mm Diameter, 70mm long, Single Use	FN-070-302	10 pcs

CE 0120

The Femcare-Nikomed Verres Needle with its red visual indicator gives tactile feedback on entry into the abdominal cavity. Safety and control at your finger tips.

- Available in two sizes: - 120mm & 150mm
- Red visual indicator
- Spring loaded stylet for tactile feedback on access to the abdomen
- Ergonomic grip
- High flow capacity

Verres Needle (Single Patient Use)



Product Description	Product Code	Pack Size
120mm Verres Needle, Single Use	FN-VN-120	25 pcs
150mm Verres Needle, Single Use	FN-VN-150	25 pcs

CE 0434

The Femcare-Nikomed Insufflation Tubing Set can be utilised during all surgical procedures where insufflation of the abdomen is required.

Insufflation Tubing Set (Single Patient Use)

Product Description	Product Code	Pack Size
 Luer Lock	FN-INS-LUER	10 pcs
 O'Ring	FN-INS-RING	10 pcs
 Click In	FN-INS-CLK	10 pcs
 15mm	FN-INS-015	10 pcs
 22mm	FN-INS-022	10 pcs
 Smoke Evacuation Tubing	FN-200-224	10 pcs

- 0.02 micron high quality filter
- 99.99% filtration efficiency
- Available with five different insufflator connectors
- Accommodates high flow rates
- Lightweight and durable



CE 0434

Femcare-Nikomed Stability Range of Trocar & Cannula Single Patient Use



Essential Abdominal Access
Made Effortless


Distributor:



FEMCARE-NIKOMED

Femcare-Nikomed Limited
Stuart Court, Spursholt Place, Salisbury Road, Romsey, Hampshire SO51 6DJ, UK
Tel: +44 (0)1794 525100 Fax: +44 (0)1794 525101
Email: enquires@femcare-nikomed.co.uk Web: www.femcare-nikomed.co.uk

Brochure Re-order Code:
BRO - 0111 (12/10) UK (N/U)



FEMCARE-NIKOMED

A FEMCARE GROUP COMPANY

Femcare-Nikomed "Stability Range" of Trocar and Cannula

The Femcare-Nikomed Single Patient Use Trocar and Cannula range is available in Shielded, Thoracic, Bladeless, Optical Bladeless and Blunt technology. The range is suitable for routine use in all minimally invasive surgery and has been designed with the surgeon and patient in mind combining performance with safety.

The Trocar and Cannula range has comprehensive features giving confidence and consistency every time:-


- Opaque cannula for easy instrument visualisation
- Easy specimen retrieval
- Silicon seal for maintaining pneumoperitoneum
- Bevelled cannula tip for easy abdominal access
- Quick & easy desufflation

Shielded Trocar & Cannula (Single Patient Use)


CE 0120

The Femcare-Nikomed Shielded Trocar and Cannula consists of a flat linear blade and a spring-loaded blade shield. The shield is designed to cover the flat linear blade once abdominal access has been gained in order to safeguard against damage to the underlying viscera.

- Reduced insertion force
- Protection from damage to the underlying viscera
- Controlled arming technology
- Spring-loaded blade shield



Available in a variety of sizes:



Product Description	Product Code	Pack Size
Shielded Trocar & Cannula with Gas Tap: 5mm diameter, 100mm long, Single Use	FN-100-105	5 pcs
Shielded Trocar & Cannula with Gas Tap: 8mm diameter, 100mm long, Single Use	FN-100-106	5 pcs
Shielded Trocar & Cannula with Gas Tap: 10mm diameter, 100mm long, Single Use	FN-100-107	5 pcs
Shielded Trocar & Cannula with Gas Tap: 12mm diameter, 100mm long, Single Use	FN-100-108	5 pcs

Blunt Trocar & Cannula (Single Patient Use)


CE 0120

The Femcare-Nikomed Blunt Trocar and Cannula allows access to the abdomen once a cut-down or Hasson technique has taken place. The blunt trocar is used to gain access to the abdominal cavity whilst the spring-loaded threaded collar is used to keep the cannula firmly in place.

- Blunt tip for decreased trauma
- For use with cut-down and open techniques
- Spring-loaded anchoring collar
- Adjustable insertion depth
- Tie down suture bars



Available in a variety of sizes:




Product Description	Product Code	Pack Size
Blunt Trocar & Cannula with Gas Tap: 10mm diameter, 100mm long, Single Use	FN-100-213	5 pcs
Blunt Trocar & Cannula with Gas Tap: 12mm diameter, 100mm long, Single Use	FN-100-214	5 pcs

Bladeless Trocar & Cannula (Single Patient Use)


CE 0120

The Femcare-Nikomed Bladeless Trocar & Cannula encompasses performance and safety. With constant downward pressure and a slight twisting motion the bladeless trocar tip separates tissue layers rather than cuts them. This results in a smaller fascial defect causing less trauma to the abdominal wall.

- Less trauma
- Smaller fascial defect
- Reduced need for post operative port site closure
- Protection from damage to the underlying viscera
- Reduced rate of port site herniation



Available in a variety of sizes:



Product Description	Product Code	Pack Size
Bladeless Trocar & Cannula with Gas Tap: 5mm diameter, 70mm long, Single Use	FN-070-215	5 pcs
Bladeless Trocar & Cannula with Gas Tap: 5mm diameter, 100mm long, Single Use	FN-100-209	5 pcs
Bladeless Trocar & Cannula with Gas Tap: 8mm diameter, 100mm long, Single Use	FN-100-210	5 pcs
Bladeless Trocar & Cannula with Gas Tap: 10mm diameter, 100mm long, Single Use	FN-100-211	5 pcs
Bladeless Trocar & Cannula with Gas Tap: 12mm diameter, 100mm long, Single Use	FN-100-212	5 pcs

Optical Bladeless Trocar & Cannula (Single Patient Use)


CE 0120

The Femcare-Nikomed Optical Bladeless Trocar and Cannula allows access to the abdomen using bladeless technology. The trocar enables new levels of control via direct vision of the abdominal layers during insertion. The bladeless tip decreases the soft tissue defect and speeds up recovery.

- Bladeless tip for reduced soft tissue defect
- Clear visualisation of abdominal layers giving more control
- Pistol grip for added comfort



Available in a variety of sizes:



Product Description	Product Code	Pack Size
Optical Bladeless Trocar & Cannula with Gas Tap: 10mm diameter, 100mm long, Single Use	FN-100-215	5 pcs
Optical Bladeless Trocar & Cannula with Gas Tap: 12mm diameter, 100mm long, Single Use	FN-100-216	5 pcs